



APPLANATION TONOMETER

DIAGNOSTIC INSTRUMENTS

MODEL **AT-1**



Aiming at new levels in quality

Applicable tonometer base



Accessories (Standard for all models)



Mounted on SM-70



Mounted on SM-30/SM-10

Specifications

Model	AT- 1 TYPE 7	AT- 1 TYPE 3	AT- 1 TYPE 4
Suitable slit lamp	SM-70	SM-30、SM-10	SM- 4 、 SM- 2
Tonometer main unit	S12-12	S12-12	S12-12
Tonometer base	S12-13	S12-14	S12-15
Measurement method	Depression type(observed with left eye on the slit lamp)		
Depression area	7.354mm ²		
Diameter of measurement depression area	3.06mm		
Measurement range	0 ~80mmHg		
Minimal scale	2 mmHg		
Measurement error	± 5 %		
Equation for conversion from mmHg to kPa	1 mmHg=101325/760Pa \div 133.3Pa 100Pa=0.1kPa		

●Kindly note that specifications may vary as improvements are carried out in the design.

TAKAGI SEIKO CO., LTD.

330-2 IWAFUNE, NAKANO-SHI, NAGANO-KEN, 383-8585, JAPAN
 TEL. +81-269-22-4511 FAX. +81-269-26-6321
 URL: <http://www.takagi-j.com> E-mail: info@takagi-j.com



CB00002 Rev.A Printed in Japan 5.2004 KY