Vantage Titan 3T is Edition

Outstanding Image quality Toshiba's unique technologies ensures high-image qualities

Toshiba implemented the most advanced Multi-phase Transmission technology that guarantees optimal Ri-distribution and homogenethy in all body regions, including addormen and cardac soans. Moreover, the new-designed Saturn Gradient coil provides enhanced image quality.



Safety Toshiba achieves true non-contrast MRA for the entire body

Toshiba's advanced non-contast MRA method allow exceptional vascular imaging without the use of contrast exception. This advanced Toshiba MR angiography suite eliminates all patient risks and provides economic benefits to your institution.



Patient comfort Large 71 on patient aperture and Planissimo" Planissimo Toshiba's advanced noise reduction to the color.

Toshiba's advanced noise reduction technology, Pantsimo, synificantly reduces the noise for all types of sequences. The 71cm of the patient aperture is the largest in the inclustry, which provides excellent comfort to patients.



Ease of Use Easy to use, Easy to learn

Our M-Power interface makes Toshiba's suite of MR software remarkably easy to learn and use, enabling all technologis to access its full range of functionalities. M-Power is a common software platform to our entire MR product line.



Eco Friendly Reducing system power consumption by "Eco Mode"

The Vantage Titan 3T utilizes a zero Helium boil-off magnet.

Moreover, the total power consumption is reduced by the "co Mode" certain electronic components go in stand-by to reduce the power consumption, but scan mode can be resumed in less than 1 second.






Vantage Titan 3T is Edition

Energy Saving Technology	Security Solution	Operation system	User Interface	Non-contrast MRA Technology	Noise Reduction Technology	Gradient Strength	Active Shield Gradient Coil	Patient Aperture	Magnet
Eco Mode ®	 Utilizing White list-based security software.⁴ Meets the requirements of Risk Management Framework (RMF), governed by the Defense Health Agency (DHA).⁸ 	Windows®7	M-Power V3.5	Time-SLIP **, FBţ FSBB, V-Trace, ASL and more advanced technologies	Plantssimo 29	30 mT/m (Standard Gradient) 45 mT/m (Saturn Gradient)	Super Slim Gradient ¹⁾	71 cm	3T High Field Homogeneity Magnet (Zero Helium Boil-off)


Vantage Titan 3T ISEASON







Advanced Clinical Solutions


The Vantage Titan 3T provides a complete suite of routine and advanced clinical solutions, including head-to-toe 2D and 3D morphological imaging with high contrast-to-noise ratio and high spatial resolution.






High-resolution MSK solution with the combination of advanced coil suites, Vantage Than 31 provides excellent high-resolution images. Moreover, advanced post-processing can be utilized.






Excellent off-center imaging
The excellent nagrets field homogeneity of Variage Titan 31 allows
uniform fits situation for off-center imaging applications such as the
shoulder or the elsow.



MPRAGE



Non-contrast CSF Dynamics Imaging

