

# Maxi-Therm®

*Single Patient-Use Hyper-Hypothermia Blanket*



**THE MAXI-THERM® BLANKET IS A DURABLE BLANKET FOR CONTROLLED WARMING AND COOLING.**

- Easy to clean, stain resistant and ideal for isolation cases
- Provides heating or cooling of the patient
- May be placed underneath or on top of the patient
- Durable materials help the blanket to last throughout the patient's entire stay

Maxi-Therm® blankets have a random flow design which results in even temperature distribution and better heating and cooling rates. The design is so efficient, the blanket can be filled with the patient in place. Set up and disposal are easy with Maxi-Therm®'s color coded, quick connectors and shut off clamps.


*Cincinnati Sub-Zero has been managing your patient's temperature needs for over 45 years.*

**CSZ**  
Cincinnati Sub-Zero

# Maxi-Therm®

*Single Patient-Use Hyper-Hypothermia Blanket*

Ordering Information		
Cat. No.	Description	Quantity
276	Adult size 24" x 60" (61cm x 152.4cm)	5/carton
274	Pediatric size 22" x 30" (55.9cm x 76.2cm)	5/carton
273	Infant size 12" x 18" (30.5cm x 45.7cm)	5/carton
Hyper-Hypothermia Units (Must purchase 286 Hoses separately to use the Maxi-Therm with any of the units listed below)		
233	Blanketrol® III	115v or 230v
222R	Blanketrol® II (not available in Europe)	230v, 240v
222S	Blanketrol® II	115v
111W	Norm-O-Temp®	100v, 230v, 240v
Reusable Connecting Hoses (must purchase separately to use blankets with hyper-hypothermia units)		
286	9' Reusable Connecting Hose (2.7m)	ea.
286-18	18' Reusable Connecting Hose (5.5m)	ea.
286-27	27' Reusable Connecting Hose (8.2m)	ea.




## Product Specifications:

### Dimensions:

Adult  
24in x 60in (61cm x 152.4cm)

Pediatric  
22in x 30in (55.9cm x 76.2cm)

Infant  
12in x 18in (30.5cm x 45.7cm)



**CSZ**  
Cincinnati Sub-Zero

Phone: 513-772-8810  
Toll Free: 800-989-7373  
Fax: 513-772-9119

[www.cszmedical.com](http://www.cszmedical.com)

© 2011 by Cincinnati Sub-Zero Products, Inc.  
99154 5/11