


NORM-O-TEMP®

Complete Hyperthermia System



NORM-O-TEMP® OFFERS MICROPROCESSOR-CONTROLLED WHOLE BODY WARMING THERAPY.

- Cost-effective heat therapy system
- Microprocessor controlled
- Pre-warming feature
- Compact size saves space

The Norm-O-Temp is a whole body system that allows conductive warming therapy to be administered. It can be used in the operating room, pre-op, recovery, or emergency department. With importance of maintaining patient core body temperature at a normothermic level during surgery becoming increasingly evident, active warming is believed to be effective in treating wounds, decreasing blood loss during procedures, and reducing post-operative infections.¹

1. Leaper, David, et al. "Maintaining Normothermia during Surgery." Hospital International. 2008. <<http://www.hospital-int.net/categories/hypothermia/maintaining-normothermia-during-surgery.asp>>

Delivering Tomorrow's Temperature
Management Solutions Today

CSZ
Cincinnati Sub-Zero

NORM-O-TEMP®

Complete Hyperthermia System

| Ordering Information | | |
|--|--|---------------------------|
| Cat. No. | Description | Quantity |
| 111W | Norm-O-Temp® | 115v, 100v, 230v, 240v |
| 118 | Stand with I.V. Pole | ea. |
| 119 | Low Profile Stand | ea. |
| Reusable Connecting Hoses (Must purchase 286 hoses separately to use blankets with Norm-O-Temp) | | |
| 286 | 9' Reusable Connecting Hose (2.7m) | ea. |
| 286-18 | 18' Reusable Connecting Hose (5.5m) | ea. |
| 286-27 | 27' Reusable Connecting Hose (8.2m) | ea. |
| Gelli-Roll® (Must purchase 286 hose separately to use Gelli-Roll with Norm-O-Temp) | | |
| 195P | Gelli-Roll® - Adult 73.5in x 21.5in x .625in (186.7cm x 54.6cm x 1.6cm) | ea. |
| 194P | Gelli-Roll® - Pediatric 31.25in x 23.75in x .625in (79.4cm x 60.3cm x 1.6cm) | ea. |
| 193P | Gelli-Roll® - Infant 20in x 13in x .625in (50.8cm x 33cm x 1.6cm) | ea. |
| Maxi-Therm® Lite Blankets (Must purchase 286 hose separately to use blankets with Norm-O-Temp) | | |
| 876 | Maxi-Therm® Lite - Adult 64in x 25in (162.6cm x 63.5cm) | 10/carton |
| 874 | Maxi-Therm® Lite - Pediatric 33in x 25in (83.8cm x 63.5cm) | 10/carton |
| 873 | Maxi-Therm® Lite - Infant 18in x 12.5in (45.7cm x 31.8cm) | 10/carton |
| PlastiPad® - reusable blankets (supplied with 9' hoses and Quick-Disconnect Couplings, except #195N & 193 CPC) | | |
| 196 | PlastiPad® - Adult 60" x 24" (152.4cm x 61cm) | ea. |
| 195N | PlastiPad® - Narrow Adult 60" x 22" (152.4cm x 55.9cm) | ea. |
| 194 | PlastiPad® - Pediatric 30" x 22" (76.2cm x 55.9cm) | ea. |
| 193 CPC | PlastiPad® - Infant (including CPC connectors) 18" x 12" (45.7cm x 30.5cm) | ea. |
| 193 | PlastiPad® - Infant 18" x 12" (45.7cm x 30.5cm) | ea. |
| 186 | 9' Blanket Extension Hose with Couplings (2.7m) | ea. |



Product Specifications:

Dimensions:

9 in x 15 in x 18.5 in
(22.9cm x 38.1cm x 46.9cm)

Weight: 33.5 lbs. (15.2 kg)

Water Temperature Range:

20°C-42°C (68°F-107.6°F)
(Heating Only)

Controller Accuracy:

Water Temp: ±.6°C (±1.08°F)

Patient Temperature:

30°C-40°C (86°F-104°F)

Safety Alarms:

Audible & Visible alarms for over and under temperature limits.

CSZ
Cincinnati Sub-Zero

Phone: 513-772-8810
Toll Free: 800-989-7373
Fax: 513-772-9119

www.cszmedical.com

CE 0344

ISO13485:2003
KEMA CERTIFIED



© 2009 by Cincinnati Sub-Zero Products, Inc.
99145 9/09